



10106924) investor in People

The Patent Office Concept House Cardiff Road Newport South Wales NP10 8QQ

PRIORITY DOCUMENT

SUBMITTED OR TRANSMITTED IN COMPLIANCE WITH RULE 17.1(a) OR (b)

REC'D 18 OCT 2000
WIPO PCT

I, the undersigned, being an officer duly authorised in accordance with Section 74(1) and (4) of the Deregulation & Contracting Out Act 1994, to sign and issue certificates on behalf of the Comptroller-General, hereby certify that annexed hereto is a true copy of the documents as originally filed in connection with the patent application identified therein.

In accordance with the Patents (Companies Re-registration) Rules 1982, if a company named in this certificate and any accompanying documents has re-registered under the Companies Act 1980 with the same name as that with which it was registered immediately before re-registration save for the substitution as, or inclusion as, the last part of the name of the words "public limited company" or their equivalents in Welsh, references to the name of the company in this certificate and any accompanying documents shall be treated as references to the name with which it is so re-registered.

In accordance with the rules, the words "public limited company" may be replaced by p.l.c., plc, P.L.C. or PLC.

Re-registration under the Companies Act does not constitute a new legal entity but merely subjects the company to certain additional company law rules.

Signed

Dated 2 October 2000

Mahoney

;25-11-99;11:46AM.; VT BY: MURGITROYD, GLASGOW eats Form 1/77 THE PAILNT OFFICE ntents Act 1977 tric 16 25 NOV 1999 The Patent Office lequest for grant of a patent FAX ice the notes on the back of this town. You can also get an Cardiff Road mlanatory leaflet from the Patent Office in help you fill in Newport Gwent NP9 1RH ds form) Your reference P230641/PKE/BOU Patent application number 25 NOV 1999 9927802.0 (The Patent Office will fill in this part) Giltech Limited Full name, address and postcode of the or of 9/12 North Harbour Estate each applicant (underline all surnames) AYR KAS SAA 000.5822001 Patents ADP number (if you know is) If the applicant is a corporate body, give the United Kingdom country/state of its incorporation "Growth Substrate" Title of the invention Murgitroyd & Company 5. Name of your agent (4) you bave one) 373 Scotland Street **GLASGOW** "Address for service" in the United Kingdom to which all correspondence should be sent G5 80A (including the postcode) 1198013 🗸 Patents ADP number (if you know tt) Date of filing Priority application number Country If you are declaring priority from one or more (day / month / year) (if you know it) earlier patent applications, give the country and the date of filing of the or of each of these earlier applications and (If you know ii) the or each application number Date of filing Number of earlier application If this application is divided or otherwise (day / month / year) derived from an earlier UK application, give the number and the filing date of the earlier application Is a statement of inventorship and of right Yes to grant of a patent required in support of this request? (Answer Yes' if: a) any applicant named in part 3 is not an inventor, or b) there is an inventor who is not named as an applicant, or c) any named applicant is a corporate body. See note (4)) Patents Form 1/77

MURGITROYD & CO. →

01633 814444;# 2

	60-11-00 11:44	01413078401	44.4049	P. 03 MURGITROYD		-2,90	Job-UZ9 01633 814444;# 3
VT F	BY: MURGITROYD, GLASGOW	;25-11-99 ; 	11:46AM ;	MUKGIIKOID	a co.	•	
			• 😀 :				<u> </u>
•	Enter the number of sheets fo following items you are filing Do not count copies of the s	with this form.	1 min man	4			
	Continuation s	ects of this form	25 NO				•
		Description	RECERSE	BY FAX			
		Claim(s)	- \	74,			
	•	Abstract	[*] (
		Drawing(*)	6 /				
0.	If you are also filing any of th	e following,					
	state how many against each		_				
	P	riority documents	•				•
	Translations of p	riority documents	·				
	Statement of inve to grant of a paten	ntorship and right					
	Request for prelim			• • •			
	and search	Patents Form 9/77	-				
	. Demiest for substi	ntive examination	· 1				
	Wedness for Ambon	(Patents Porm 10/77					
	Ąn	y other document Olean specifi					
11.			I/We re				easis of this application
			Signatu	ie meyos	,		Date 25.11.1999
				MURGITROYD	& COMPA	NY·	
12	Name and daytime telephor person to contact in the Un	ne number of ited Kingdom	Ве	verley Ouzma	n . 01	41 307	8400
Af	terning ter an application for a pain communication of the investil be informed if it is necess nited Kingdom, Section 23 of ritten permission from the I inited Kingdom for a patent ommunication has been give	ntion spoura be p ary to probibit or f the Patents Act Patent Office unles for the same inve	restrict your in 1977 stops you is an application intion and elibe	ivention in this from applying f in has been filed or no direction f	way. Purt or a patei I at least b	bermore it abroa weeks	r, if you live in the A without first getting beforeband in the
N	otes) If you need belp to fill in	this form or you i	ave any quest	ons, please cont	act the Pa	sent Off	ice on 0645 500 5 05.
h) Write vour answers in ca	ottal letters using	black ink or yo	u may type tber	773.		
c)) If there is not enough spa sheet of paper and write attached to this form.	Sam all also males	ant datable on	any part of this	form, blee	ase conti nuation	inue on a separate sbeet sbould be
d) If you bave answered Ye						
8							
f) For details of the fee and	ways to pay plea	pe contact the E	atent Office.			
	•	ı					Patents Forma 1/

Job-029

R-290

P.03

1

GROWTH SUBSTRATE

10 11

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

The present invention provides a growth substrate for cell culture. In particular, the present invention provides a cell culture growth substrate for tissue engineering.

Tissue engineering is expected to transform orthopaedics treatments, cancer therapy and the treatment of chronic degenerative diseases. engineering concerns the provision of a graft comprising living cells or suitable substrate to sustain the growth of such cells which integrate into the patient providing expedited wound healing and repair or an alternative drug delivery or gene therapy delivery system. The tissue engineering graft may be an autograft, allograft or xenograft. Autografts are formed with the patient's own cells, cultured with a suitable growth medium or substrate. Allografts rely upon cells donated from an alternative same species source (including cadaver or foetal sources) whilst xenografits rely upon cells donated from other species. Both allografts and xenografts may be treated to

32

01633 814444;# 6

2

minimise autoimmune rejection of the graft following 1 implantation. 2 3 There are numerous potential applications for tissue 4 engineered grafts, including reconstructive surgery, 5 orthopaedics or dental applications, burn treatments 6 or ulcer treatments (including venous ulcers and 7 diabetic foot ulcers). A number of tissue engineered В grafts have been described in the literature (see 9 Dutton, "Tissue Engineering", Genetic Engineering 10 News, Vol 18, No 8, April 15, 1998). 11 12 Examples of such tissue engineered grafts include 13 APLIGRAF (Trade Mark) which is a bilayer graft. 14 including both differentiated keratinocytes and a 15 layer of fibroblasts in a collagen matrix. APLIGRAF 16 has been used as a skin graft, particularly for burns, 17 diabetic foot ulcers, excisional surgery and venous 18 ulcers (Bender, *Healing of Difficult to Heal Wounds 19 Using a Bilayered Skin Construct", 11th Annual 20 Symposium on Critical Issues in Surgery-Wound Healing, 21 Science and Technology, 3-5 December 1998, St Thomas, 22 23 US Virgin Islands). Other bioengineering skin equivalents include INTEGRA (Trade Mark), a xenograft 24 of-bovine collagen, glycosaminoglycans (GAG) and 25--silastid sheet; ALLODERM (Trade Mark), an allograft of 26 treated cadaver skin; and DERMOGRAFT (Trade Mark), an 27 allografit of neonatal fibroblasts on a polyglactin 28 Tissue engineered grafts for bone include 29 scaffold. RAINBOW (Trade Mark) of IsoTis BV which is a 30 biomimetic coating which allows a bone-like layer to 31

grow over metal prosthesis and serves as a scaffold

for bone growth, and also EMBARC (Trade Mark) which is a resorbable bone repair material.

3

4 Despite the numerous tissue engineering grafts

5 currently being developed, there is still a demand for

6 further and improved products. We have now found that

7 water-soluble glass acts as a support or matrix for

8 cell growth and hence the glass has utility in tissue

9 engineering.

10

13

11 The present invention thus provides a cell culture

12 growth substrate comprising a water-soluble glass

matrix adapted to sustain the growth of living cells.

14 Preferably the substrate will comprise or have at

15 least a portion of the surface thereof coated with

16 living cells.

17

18 In one embodiment the cell culture growth substrate is

19 pre-seeded with living cells and hence the matrix

20 comprises or has at least a portion of its surface

21 coated with living cells.

22

23 In one embodiment, the cell culture growth substrate

24 will be useful as a tissue graft, i.e. is designed for

25 implantation into a patient to replace or promote

26 repair of damaged tissues.

27

28 The water-soluble glass matrix will of course be

29 biocompatible. Generally, the biodegradation of the

30 water-soluble glass following implantation of the

31 graft into a patient will be pre-determined to be

25-11-99 11:44 01413078401 P.08 R-290 Job-029
3ENT BY:MURGITROYD, GLASGOW ;25-11-99 ;11:47AM ; MURGITROYD & CO.→ 01633 814444;# 8

4

compatible with the timescale required for regrowth of the tissues concerned.

3

12 .

The glass present in the graft acts as a cell support
matrix and will function as such in vivo. Thus the
graft can be used directly in vivo to provide a
temporary biodegradable scaffold which will encourage
ingrowth of surrounding tissues. In other embodiments
pre-seeding of the graft with a pre-selected cell
type, and optionally growth of that cell type, prior

11 to implantation may be desirable.

In an alternative embodiment, the cell culture growth 13 substrate is intended for non-clinical purposes, for 14 example in bio-reactor and fermentation technologies 15 for the production of drugs and other biologically 16 Organisms usually grow with derived chemicals. 17 increased confluence on surfaces, and enzyme reactions 18 (and many other biochemical reactions) are generally 19 most efficient when the enzyme is bound to a reaction 20 Beads, sinters and fibres can be used to 21 surface. provide the required mechanical support, with large 22 (productive) surface areas and additional features 23 such as controlled inorganic micro-nutrient supply, 24 contamination control, pH buffering and a 25 biocompatible carrier which will allow the subsequent 26 transfer or filtration of cells, enzymes or other 27 components bound to its surface on completion of the 28 reaction stage. 29

30

Conveniently the water-soluble glass matrix may be in the form of water-soluble glass fibres and reference BY: MURGITROYD, GLASGOW

30

31

MURGITROYD & CO.→

is made to our WO-A-98/54104 which describes the 1 production of suitable glass fibres. Whilst the glass 2 fibres can be used in the form of individual strands, 3 woven (e.g. a 1 x 1 basket weave) or non-woven mats 4 may also be produced from the fibres and used as the 5 The individual fibres of a non-woven mat may 6 be gently sintered together to obtain coherence of the 7 Alternatively, the fibres may be used as strands. 8 glass wool and this form of matrix is especially 9 suitable where the graft requires a 3D shape. 10 11 Alternatively, the water-soluble glass matrix may be 12 produced from finely comminuted glass particles (for 13 example having an average diameter of from 50 µm to 6 14 Optionally, the glass particles may be sintered 15 together to form a porous structure into or onto which 16 cells may be seeded and in this embodiment the glass 17 particles will have a preferred diameter of from 53 µm 18 to 2 mm, preferably 400 μm to 2 mm. Again, a three-19 dimensionally shaped graft may be produced (if 20 necessary individually tailored to be compatible with 21 the wound site of the patient) from the sinter. 22 Alternatively, particles following a Fuller curve 23 packing distribution and having a range of diameters 24 of 0.3 mm to 5.6 mm may be used. 25 26 In a further embodiment the glass may simply be in the 27 form of a glass sheet, which may be substantially 28 planar or may be contoured to a required shape. 29 Etched, ground or patterned glass sheet may be used in

addition to plain surfaced glass.

25-11-99 11:44 01413078401 P.10 R-290 Job-029 ENT BY:MURGITROYD, GLASGOW ;25-11-99 ;11:48AM ; MURGITROYD & CO.→ 01639 814444;#10

6

The water-soluble glass preferably includes 1 phosphorous pentoxide (P_2O_5) as the glass former. 2 3 Generally the mole percentage of phosphorous pentoxide 4 in the glass composition is less than 85%, preferably 5 less than 60% and especially between 30-60%. 6 7 Alkali metals, alkaline earth metals and lanthanoid 8 oxides or carbonates are preferably used as glass 9 modifiers. 10 11 Generally, the mole percentage of alkali metals, 12 alkaling earth metals and lanthanoid oxides or 13 carbonates is less than 60%, preferably between 40-14 15 60%. 16 Boron containing compounds (e.g. B203) are preferably **17**. used as glass additives. 18 19 Generally, the mole percentage of boron containing 20 compounds is less than 15% or less, preferably less 21 than 5% 22 23 Other compounds may also be added to the glass to 24 modify its properties, for example SiO2, Al2O3, SO3, 25 sulphate ions (SO42-) or transition metal compounds 26 (e.g. first row transition metal compounds). 27 28 Typically the soluble glasses used in this invention 29 comprise phosphorus pentoxide (P2O5) as the principal 30 glass-former, together with any one or more 31

glass-modifying non-toxic materials such as sodium

30

31

32

oxide (Na_2O) , potassium oxide (K_2O) , magnesium oxide 1 (MgO), zinc oxide (ZnO) and calcium oxide (CaO). 2 rate at which the glass dissolves in fluids is 3 determined by the glass composition, generally by the 4 ratio of glass-modifier to glass-former and by the 5 relative proportions of the glass-modifiers in the 6 By suitable adjustment of the glass 7 glass. composition, the dissolution rates in water at 38°C 8 ranging from substantially zero to 25mg/cm²/hour or 9 more can be designed. However, the most desirable 10 dissolution rate R of the glass is between 0.001 and 11 2.0mg/cm²/hour. 12 13 The water-soluble glass is preferably a phosphate 14 glass, and preferably comprises a source of metal ions 15 which confer either antimicrobial protection or 16 enhanced cell growth, or both, or which are useful 17 trace elements. Examples include silver, copper, 18 magnesium, zinc, iron, cobalt, molybdenum, chromium, 19 manganese, cerium, selenium, and these metal ions can 20 be included singly or in any combination with each 21 Where silver ions are of interest, these may other. 22 advantageously be introduced during manufacture as 23 silver orthophosphate (Ag3PO4). The glass preferably 24 enables controlled release of metal ions and other 25 constituents in the glass and the content of these 26 additives can vary in accordance with conditions of 27 use and desired rates of release, the content of 28 silver generally being up to 5 mole %. While we are 29 following convention in describing the composition of

the glass in terms of the mole % of oxides, of halides

and of sulphate ions, this is not intended to imply

25-11-33 11:44 01413078401 P.12 R-290 Job-029

ENT BY: MURGITROYD, GLASGOW ;25-11-99 ;11:48AM ; MURGITROYD & CO. → 01633 814444;#12

8

that such chemical species are present in the glass 1 nor that they are used for the batch for the 2 preparation of the glass. 3 4 The optimum rate of release of metal ions into an 5 aqueous environment may be selected by circumstances 6 and particularly by the specific function of the 7 released metal ions. The invention provides a means 8 of delivering metal ions to an aqueous medium at a 9 rate which will maintain a concentration of metal ions 10 in said aqueous medium of not less than 0.01 parts per 11 million and not greater than 10 parts per million. .12 . some cases, the required rate of release may be such 13 that all of the metal added to the system is released 14 in a short period of hours or days and in other 15 applications it may be that the total metal be 16 released slowly at a substantially uniform rate over a 17 period extending to months or even years. 18 particular cases there may be additional requirements, 19 for example it may be desirable that no residue 20 remains after the source of the metal ions is 21 exhausted or, in other cases, where the metal is made 22 available it will be desirable that any materials, 23 other than the metal ions itself, which are 24 simultaneously released should be physiologically 25 harmless. In yet other cases, it may be necessary to 26 ensure that the pH of the resulting solution does not 27 fall outside defined limits. 28 29 Generally, the mole percentage of these additives in 30 the glass is less than 25%, preferably less than 10%. 31

;25-11-99 ;11:49AM ;

9

01633 814444;#13

- The cells may be any suitable cells required for 1
- Particular mention may be made of 2 grafts.
- keratinocytes, fibroblasts, chrondrocytes and the like 3
- as preferred cell types. Mention may also be made of 4
- stem cells (mesenchymal, haematopoetic, and 5
- embryonic), Schwaan cells, keratinocytes (epithelial 6
- cells), chondrocytes, osteoblasts, endothelial cells 7
- and other fibroblasts, cardiac cells (and other 8
- myoblasts), pancreatic β cells and periodontal tissues 9
- e.g. Dentine, but the invention is not limited to 10
- these cell types alone. 11

12

- Embodiments of the invention will be described with 13
- reference to the following non-limiting examples and 14
- Figures in which: 15

16

- Chondrocytes forming a monolayer on a glass Fig. 1 17
- fibre (Example 1) as viewed by laser scanning 18
- 19 confocal microscope.

20

- Fluorescent microscopy of HUE cells on MATT01 21 Pig. 2
- glass fibres (see Example 2). 22

23

- Fluorescent microscopy of HUE cells on MATT04 24 Fig. 3
- glass fibres (see Example 2). 25

26

- SEM picture of L929 cells on glass surface at 27 Fig. 4
- x30 magnification (see Example 3). 28

29

- SEM picture of L929 cells on glass surface at 30 Fig. 5
- x170 magnification (see Example 3). 31

25-11-99 11:44 01413078401 P.14 R-290 Job-029 ENT BY:MURGITROYD, GLASGOW ;25-11-99 ;11:49AM ; MURGITROYD & CO.→ 01633 814444;#14

10

Fig. 6 SEM picture of L929 cells on glass surface at x215 magnification (see Example 3).

Fig. 7 SEM picture of L929 cells on glass surface at x610 magnification (see Example 3).

NT RY:MURGITROYD, GLASGOW ;25-11-99;11:49AM; MURGITROYD & CO.→ 01633 814444;#15

11

1 Example 1

2

Introduction

4

5 Controlled Release Glass (CRG) is a phosphate-based

6 material which degrades at a predeterminable rate.

7 The potential for using CRG as a cartilage engineering

8 matrix has been assessed using isolated equine

9 chondrocytes with in-vitro techniques. The glass was

10 provided in fibrous form in three different

11 compositions. The three CRG compositions provided

12 have showed potential as a tissue engineering

13 substrate.

14 15

Materials and Method

16

27

17 A total of 200,000 chondrocytes isolated from horse

18 articular cartilage were added to each 2 cm well in a

19 24 well plate. Every well contained 0.02 grams of

20 glass fibre sample. Four different fibres F1 to F4

21 (diameters 20-30 μm) were analysed: F1 - containing

22 Fe₂O₃ and NaF, F2 - containing Ce₂O₃ and Se. The

23 composition of glasses used to form F1 to F4 are set

24 out below in Table 1. The culture medium (containing

25 10% FCS was changed daily. At time periods of 3

26 days, 1 week and 2 weeks, the samples were stained

using rhodamine phalloidin and oregan green for the

viewing of actin and tubulin using a laser scanning

29 confocal microscope. At the same time periods, the

30 cell supernatant was removed and stored at -80°C until

31 analysis on cell viability and type II collagen

32 production could be performed. Production of type II

25-11-99 11:44 01413078401 P.16 R-290 Job-029
ENT BY:MURGITROYD, GLASGOW ;25-11-99 ;11:49AM ; MURGITROYD & CO. → 01633 814444;#16

12

collagen was analysed by using RT-PCR analysis on the 1 cDNA from the chondrocyte population in contact with 2 the glass fibres. The total RNA was prepared from the 3 cell population by the addition of 1 ml of TRIzol 4 (SIGMA) to the cell population for 5 minutes. After 5 this time, the TRIzol was retrieved and stored at -6 80°C until RT-PCR analysis could be carried out. 7 RT-PCR analysis was performed by tagging with primers 8 for collagen type II and with gapDH for cell 9 10 viability.

11

2ymography was also performed at time periods of 4
days, 1 week and 2 weeks for detection of matrix
metalloproteinases (MMP's) produced by the
chondrocytes.

		BELEC
н		RECORD
Table		BATCH
~	C	œ

1				,	_	T	7	
Physical	Form		F,C + R	F,C + R	F.R + C		× ·	
Solution Ratus		Non-Annealed @ 37.5°C (mg.cm ⁻¹ .hx ⁻¹)	N/A	N/A	W/W		0.5614	
Solution		TOTAL Americal @ 17.5°C (mg.cm ⁻² .hr ⁻¹)	N/A	N/A	N/A		0.351	
		TOTAL	9 <u>1</u>	99.6		ı	BI	
					26 6 30 1 7 1			
		Çelo		2.6	3			
		Map	1:1		5	:		
	formulation as moles	Pejo, Haf Celo,	1.76	56.0		26.5		
Boton		1	5.75 1.5	5		1:43		
tion as		MgO KiO BiO: MuO	5, 75	2		2.00	2	
I I		N ₂ O					ន	
		Sign Sign					4	
		P ₂ O ₅	4K 78		100	45.05	46	
		yeio					3	
		CaO	19 96	29.74	17.94	17.03		
		Na ₃ O	20 20	65.53	25.7B	25.19	32	
	}		į		2.2	E.	F4	L

F=fibres; C=cullet; R=rods 9 1

01413078401 Job-029 25-11-99 11:44 P.18 R-290 01633 814444;#18

14

ENT BY: MURGITROYD, GLASGOW ;25-11-99 ;11:50AM ; MURGITROYD & CO. →

> Regults and Conclusions 1

3 Chondrocytes adhered to all three types of fibre

At the 3 day time period, the cells appeared 4

5 to be rounded. At 1 week and 2 weeks, confocal

microscopy indicated cell proliferation between all 6

7 time periods. At 1 week and 2 weeks, the cells were

8 elongated and formed a monolayer along the fibre

length as can be seen in Figure 1. 9

10 The RT-PCR analysis showed that fibres F2 and F3 were 11

12 . producing collagen type II.up to and including the two

13 week time period indicating that the cells retained

14 their chondrocytic phenotype.

15

2

The zymography performed on F2 and F3 showed that the 16

17 cells in contact with these fibres produced MMP2 at

18 all three time periods, but in a greater quantity at 2

19 weeks than 1 week, and at 1 week than 4 days.

20 increase of MMP2 production is expected, as the cells

21 were seen to have increased in number at these time

22 periods from the confocal microscope analysis.

23

24 In conclusion, all three fibres types showed cell

25 adherence and the chondrocytes adhered to F2 and F3

appear to retain the ability to produce type II 26

27 collagen.

Example 2 1 2 Biological Evaluation of Non-woven Mat Fibres 3 4 1. Objective 5 6 Using in-witro techniques determine: 7 The cytotoxicity of a series of five nona. 8 woven mat CRG fibres. 9 The potential of the fibres as a cell b. 10 substrate matrix. 11 12 13 2. Scope 14 The test procedures apply to all fibre samples. 15 16 3. Equipment and Materials 17 18 3.1 Equipment 19 Laminar air flow hood 20 3.1.1 Incubator maintained at 37°C/5% carbon 3.1.2 21 dioxide 22 Refrigerator at 4°C 3.1.3 23 3.1.4 Freezer at -18°C 24 3.1.5 Vacuum source 25 Phase contrast microscope 3.1.6 26 27 3.2 Materials 28 Sterile plastic-ware pipettes 3.2.1 29

Sterile glass pipettes

24 well Sterile dishes

30

31

3.2.2

3.2.3

1	3.2.4	Surgical grade forceps
2	3.2.5	Surgical grade scissors
3	3.2.6	Sterile Universal containers
4	3.2.7	L929 cell culture line (ATCC NCTC Clone 929)
5	3.2.8	Human Umbilical endothelial cells (primary
6		cell source, Liverpool Women's Hospital)
7	3.2.9	TCPS negative control
8	3.2.10	CRG fibres:
9		D021298F1 (MATT01)
10		D301198F1 (MATT02)
11		D100299F1 (MATT03)
. 12		D161298F2 (MATT04)
13		D171298F2 (MATT05)
14		All CRG fibres were supplied non-sterile in
15		quantities 8g-38g. The compositions of CRG
16		fibres used (MATT01 to 05) are set out below
17		in Table 2.

Table 2

BATCH RECORD SELECTION

Physical Form # 37.5°C (mg.cm⁻³.hr⁻³) Non-America N/A N/A 0.0151 0.0165 0.02 Solution Bates Annealed @ 37.5°C (mg.cm⁻³.hr⁻³) N/A N/A 0.0095 0.0143 0.0177 100 100 100.47 100 2.27 1.96 CejO 1.54 0.96 76,05 Formittion as molet 3.68 6.36 K₁O B₂O₃ B 8 8 50 50 47.04 50.2 O'GY 25.19 17.03 27.98 8 8 MATTO1 MATTO2 MATTO4 MATTO6 900

Rerods; Fefibres; Cecullet

1	4.	Proc	edure	
2				
3		4.1	Test samp	ole preparation
4			4.1.1	Test samples were cut to the
5				appropriate size (see section
6				4.2.1).
7			4.1.2	Tissue culture polystyrene was
8				employed as a negative control.
9				The controls were not in the same
10				physical form as the test
11				material.
12				
13		4.2	Fibres we	ere examined in contact with the
14		L929	cell line	before any cleaning procedure.
15		Fibre	s were ex	samined in contact with both cell
16		lines	after cl	eaning in acetone, washing in PBS
17		and s	 sterilisin	ng in a dry oven at 190°C for 2
18		hours	 \$.	
19				
20		4.3	Cell prep	paration
21			4.3.1	A cell subculture was prepared 24
22				hours before being introduced to
23				the fibres.
24				
25		4.4	Test proc	edure
26			4.4.1	A small "bed" of the fibres was
27				placed in the bottom of each well.
28			4.4.2	The cell/medium preparation was
29				gently pipetted onto the fibre
30				bed.

	•	
1	4.4.3	The 24-well plates were incubated
2		and examined at 24 hours and 48
3		hours.
4		
5 4.5	5 Interpre	tation of results
6	4.5.1	At the conclusion of the
7		incubation period the plates are
8		removed from the incubator and
9		examined under phase contract
10		microscopy using x10 and x20
11		objective lenses.
12	4.5.2	Each test and control material was
13		initially evaluated using the
14		scoring system detailed below.
15		This evaluation was based on the
16		appearance of the cells which were
17		attached to the TCPS surface. It
18		was not possible to carry out such
19		an evaluation on the cells
20		adhering to the fibres.
2.7	l	

Table	3	:	Reactivity	Responses	

Grade	Reactivity	Conditions of all cultures
0	None	Discrete intracytoplasmic granules; no cell lysis
1	Slight	Not more than 20% of cells are round, loosely attached and without intracytoplasmic granules; occasional lysed cells are present
2	Mild	Not more than 50% of the cells are round and devoid of intracytoplasmic granules; extensive cell lysis and empty areas between cells
3	Moderate	No more than 70% of the cell layers contain rounded cells and/or are lysed
4	Severe	Nearly complete destruction of the cell layers

25-11-99 11:44 01413078401 P.24 R-290 Job-029

01633 814444;#24

SENT BY: MURGITROYD, GLASGOW ;25-11-99 ;11:51AM ; MURGITROYD & CO.→

4.6 Cytotoxicity Results

The following table highlights the results obtained following two separate tests. Two or four readings were taken at each test. In all cases negative control (TCPS) provided a 0 grade.

8 Table 4

Material Code	Grade Test 1 L9292	Waterial Code	Grade Test 2 1929	Test 2 HUE
MATT01	0	MATT01	0	0
MATT02		MATT02	0	0
MATT03	_	MATT03	-	0
MATT04	0	MATT04	0	0
MATT05	0	MATT05	0	0

Comments

The results as detailed provide a very subjective assessment of material cytotoxicity. Where a grade 0 is shown, there was no evidence of toxicity and a confluent healthy monolayer of cells was present. Where there was evidence of contamination or where the cell monolayer is difficult to evaluate no score has been given.

4.7 Cell Substrate Results

The following table (Table 5) details the cellfibre interactions and general cell culture
conditions observed by phase contrast microscopy.

As stated before phase-contrast images of the cells
on the fibres are poor. A staining procedure was

carried out with the HUE cells. This procedure 1 uses a #luorescent staining technique (ethidium 2 bromide and acridene orange) to identify cell 3 viability. All observations were after 48 hour 4 contact between cells and fibres. 5

6 7

Table 5

		<u> </u>	120 TO 120 TO 1	MATT03	MATTO4	MATTO5
	MA Cells viabl	1000		Contamination	viable. Healthy monolayer. There is some evidence of cell adhesion onto the fibres.	Cells are viable. Healthy monolayer. There is some evidence of cell adhesion onto the fibres. Cells viable.
1929 sterile fibres		pH levels are low. Cells are viable. There is no obvious cell adherence	pH is low. There seems to be evidence of contamination- though this may be degrading glass. Difficult to make any comment on cell viability.	monolayer. There is some evidence of coll adhesion	Healthy monolayer. There is some evidence of cell adhesion onto the fibres.	
HUR sterile fibres	viab thous gran appe The pH h drop Some can adhe	le gh ular in arance. medium	cells are viable. There is no obvious cell adherence to the fibre. contami though be degr glass. Difficu	pH is low. There seems to be evidence of contamination- though this may be degrading	the control wells. There is some evidence of cell attachment	difficult to observe by phase

The images below were obtained following the vital 9 staining procedure and examined by fluorescent 10

microscopy. 11

12

13

14

15

16

17

As well as demonstrating cell viability the procedure permitted a better evaluation of the cells attaching to the fibres. The cell-fibre interaction was much better than that indicated by phase contrast microscopy. It was noted that

P. 26 R-290 Job-029 01633 814444;#26

01413078401 ENT BY: MURGITROYD, GLASGOW

25-11-99 11:44

32

;25-11-99 ;11:51AM ; MURGITROYD & CO.→

22

MATT04 and MATT05 had excellent cell adherence. 1 MATT01 permitted a good cell adherence. There was 2 cell attachment with MATT02 and MATT03 although 3 this was poor in comparison with 01, 04 and 05. 4 5 Example 3 6 7 A cell suspension (in complete cell culture medium 8 supplemented with 5% foetal calf serum) at a 9 concentration of approx. 5 x 105 cell/ml was 10 introduced to an established mouse fibroblast cell 11 line (L929). 12 13 The material/cell interaction was examined using 14 phase contrast microscopy at 24, 48 and 72 hours. 15 In particular the following materials were examined 16 (see Table 6 for composition of the glasses 17 referred to by batch number). 18 19 a) Glass sheet (flat); code 1051098-1 20 21 Cells can be seen adhering to the material and 22 remain in contact with the material following 23 The cell sequential transfer between dishes. 24 morphology is rounded and the growth rate is 25 considerably slower than observed with cells on the 26 Nevertheless there is evidence of control dishes. 27 cell division taking place on the surface. 28 29 30 31

b) Sintered glass beads (smooth surface); code BX
D221098-1, Sintered glass beads (rough surface);

code BX-D221098-1

4

It is more difficult to make the observations with these samples using phase contrast. However, cells are clearly present on the surface of both rough and smooth samples. The cell population is certainly increasing with time up to the 72 hour period. Again, this is following sequential

11 transfer at 24 hours.

Table 6

Batch	For	molet		TOTAL	Solutio	Physical Form	
Number	Na ₂ O	CaO	P205	1	Annealed e37.5°c (mg.cm'2.hr'1	Annealed e37.5°c (mg.cm'2.hr'1	
1051098-1	25	28	47	100	0.0991	0.1364	R+5
D221098-1	11	42	47	100	0.0377	0.0446	G+R

G=GRANULES R=RODS S=SHEETS

Sample SEMs were obtained (see Figs. 4 to 7) after cells had been in contact with the glass for 72 hours, fixed in 2.5% glutaraldehyde and dehydrated with alcohol. The samples were gold coated before viewing. The magnification is indicated on Figs. 4 to 7.

T BY: MURGITROYD, GLASGOW, ;25-11-99 ;11:52AM; MURGITROYD & CO.→

01633 814444;#28

Figure 1




Figure 2




FIGURE 1 [MATTO1]

The bright areas represent viable cells [HUE]. The image shows an area with bundles of fibres radiating in many directions. In most cases the cells are rounded and not elongated on the fibres.

Figure 3




FIGURE 2 [MATT04]

The bright areas represent viable cells [HUE]. Cells can be seen elongated on the fibres. In this image most of the fibres are oriented in the same direction. There is excellent cell coverage. This image is also representative of the result obtained with MATT05

Comment

Of the five fibre compositions examined MATT04 and MATT05 are providing an excellent substrate for cell adhesion. MATT01 has large numbers of cells adhering although the cell morphology is more rounded than that seen on the control surface. MATT02 and MATT03 show cells adhering but in much reduced numbers. There is no evidence of cytotoxicity with any of the fibres examined.

Figure 4




Figure 5




Figure 6





Figure 7



Pet NO: 6600 / 03424

Phil Tream

GB0003424

This Page is Inserted by IFW Indexing and Scanning Operations and is not part of the Official Record

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

D	efects in the images include but are not limited to the items checked:
	☐ BLACK BORDERS
	☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
	☐ FADED TEXT OR DRAWING
	☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
	☐ SKEWED/SLANTED IMAGES
	☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
	☐ GRAY SCALE DOCUMENTS
	☐ LINES OR MARKS ON ORIGINAL DOCUMENT
	☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
	·

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.